

Thermo Evolution 220 UV-Vis

Sample Prep

1. Solution: Prepare a 1×10^{-4} M solution in CHCl_3
2. Thin Film: Spin-cast on glass slides from 2 mg/mL solutions in chlorobenzene

Running a Sample

1. Turn on the UV-Vis (switch is on the back). The lamp doesn't need to warm up :)
2. Turn on the computer.
3. Open the INSIGHT 2 program and select the method you want (e.g. scan)
4. Click "Settings" and set the desired parameters:
 - a. **Measurement tab**
 - i. Method template description
 - ii. Baseline correction (if necessary)
 - b. **Instrument tab**
 - i. Data mode (absorbance or transmittance)
 - ii. Wavelength range for measurement
 - iii. Integration time, interval, scan speed, etc.
 - c. **Sample tab**
 - i. Number of samples and base name for files
 - d. **Peak Pick tab**
 - i. Sensitivity for automatic data analysis (can also do this manually later)
5. Go to "file (menu) > save workbook settings as template" to create a template
6. Click "measure"
7. Insert blank--make sure to close the cover
8. click "OK" to measure background
9. "Confirm sample list" box will pop up; verify the sample information and insert the first sample in the batch

10. click “continue” to take measurements
11. Pick peak(s) if you didn’t choose to auto-pick them by going to “analyze > peak pick”
12. Click “reports” to save or export the data
 - a. Select sample data from “samples” tab
 - b. Select the order top include the data in the “reportable data” tab
 - c. Use the “layout” tab to pick which items to display
 - d. Preview the report if necessary
 - e. Export by clicking “export”
13. When you're finished, close the program
14. Turn off the computer and UV-Vis